

# PowerGrip Bipolar Laparoscopic Instruments





## PRODUCT HIGHLIGHTS

- PowerGrip bipolar laparoscopic instruments facilitate precise dissecting, grasping, cutting and coagulation
- The PowerGrip bipolar handles and inserts are reusable
- Jaw mechanism enables adjustable opening and closing of the jaws with very high pressure while grasping and cutting
- Smooth handle operation provides for an optimal tactile feel and mechanical advantage
- Instrument jaws are insulated to the end of the grasping and cutting zone in order to avoid unintentional coagulation
- Rotation knob allows electrode insert to be rotated up to 360° with the forefinger
- Modular instruments facilitate ease of assembly and cleaning

## REUSABLE POWERGRIP HANDLE ASSEMBLY

		Catalog #		
Handle Assembly			34cm	45cm
	PowerGrip Handle, 2-Pin	BP7500	---	---
	PowerGrip Shaft	---	BP7550	BP7550-45

## REUSABLE POWERGRIP BIPOLAR INSERTS

		Catalog #		
Tip Style			34cm	45cm
		Bipolar Forceps	BP7601	BP7601-45